

OXYGEN CONCENTRATOR

MODEL: DMZ-10



Single flow



Dual flow



USA imported humidifier
6PSI pressure release



Innovative "cooling"
technology to improve
reliability and life of the unit



Big LCD display
24H*365D work continuously

Specification

Flow Rate	1-10LPM
Purity	93% ($\pm 3\%$)
Outlet Pressure	0.04-0.07MPa(6-10PSI)
Sound Level	≤ 50 db
Power Consumption	≤ 880 w
LCD Display	Switch times, operating pressure, present working time, accumulated working time, presetting time from 10mins to 40hours
Alarm	Power Failure alarm, High&low pressure alarm
Gross Weight	29 Kgs
Package Size	435×435×640mm